

# Catalyst Testing in FCC Naphtha Hydrotreating

Key insights for smart selection

Cansu Mai | hte  
Dan Miskin | Evonik

# Who we are: hte in numbers

hte is the world leading solution provider for lab-scale R&D workflows



Clients from **38 countries** with **> 250 systems** delivered all over the world



Largest high throughput catalysis lab worldwide with **> 50 reactor units** and up to **1,000 reactors** running 24/7



**> 350 employees** with highly skilled scientists and engineers



**Founded in 1999** in Heidelberg, Germany and **25 years of experience**




Financially sound & reliable ownership structure with **BASF since 2008**

# Our approach

---



Downscaling  
commercial  
processes



Parallelization  
&  
automatization



Digitalization



High quality &  
commercially  
relevant data



# Selected customer portfolio

## Oil & Gas



## Refineries



## Chemical




## Catalyst



## Battery & Electrolysis



# hte technology for catalyst & process development



# 16/24-fold trickle bed reactor unit



## Flexible tool for testing hydrocracking catalysts



<b>Area of application</b>	HDS, HDN, HDC, HDO, HDA, HDM, Hydrosom., Dewaxing, lubes, bio-feeds	<b>Parallelization</b>	16 individual reactors or 8 x 2 reactors in series (= 8 x 3 reactors); interstage sampling and dosing
		<b>Temperature</b>	up to 450°C 16 individual heating blocks
		<b>Pressure</b>	15 to 260 bar
		<b>Catalyst volume</b>	1-10 mL (Extrudates)
		<b>Feed requirements</b>	~200 L/month operating 16 parallel trains
		<b>Product collection</b>	~1 L of TLP for each reactor/week
		<b>Feeds</b>	H2 feed gas, N2 make-up gas 1 or 2 liquid feeds Diesel, VGO, resid, bio-based feed, etc.
		<b>Operation mode</b>	Down or up flow, once-through, in series
		<b>Online analysis</b>	Gas phase HC, H2, H2S, CO, CO2
		<b>Offline analysis</b>	SIMDIST, density, S, N...

# Trickle bed test unit & reactor loading


## General set-up



- Catalyst extrudates of different sizes tested in:
  - 15x 4.80mm and 1x 6.35mm trickle bed reactors
  - Individual capillary heating to adjust the pressure drop
- 12 Naphtha hydrotreating catalysts and 3 Selective hydrogenation catalysts evaluated in parallel
  - Individual reactor heating system
  - Reactor temperature differences up to 150°C

# Project timeline

- **6-8 weeks prior**
  - Alignment on test protocol and reactor loading
  - Catalyst and feed shipments
- **After stage**
  - Final report issued
  - Shipment of samples of interest





# Case study - FCC naphtha hydrotreating

## Balancing HDS activity with octane retention

- High HDS severity risks excessive hydrogenation, leading to octane loss and increased H<sub>2</sub> demand.

- Project specific advanced analytics

✓ Workflow for reactive sulfur species analysis: selective oxidation prep and SCD-based speciation

### ✓ Octane retention screening

Bromine number = *total unsaturation*

Maleic anhydride titration = selective reactivity of octane-contributing **olefins** and **dienes**

**Sulfur speciation  
&  
Octane number**

**HDS vs. Olefin  
saturation**

**Different  
challenges**

**Feed stability  
& Clogging:  
Dienes, olefins,  
thiols**

✓ Operating conditions tuned to moderate temperatures (~140–310 °C) and optimized H<sub>2</sub>/oil ratios (~320 Nm<sup>3</sup>/m<sup>3</sup>.h) to minimize hydrogen use.

✓ Catalyst selection focused on high selectivity for thiophenic sulfur removal while preserving *olefins*.


- Reactive feed handling

✓ Condenser temperatures were kept low (40–50 °C) to stabilize vapor–liquid equilibrium and minimize light-end losses

# Results


## Mass balance & RON

Mass balance vs. TOS



- Carbon mass balance accounted across gas and liquid phases, resulting in 100% ( $\pm 3\%$ ) for most datapoints collected

RON vs. TOS




- RON calculated from GC composition analysis
- Bromine number and maleic anhydride titration also performed throughout experiment conditions

# Results


## Hydrogen consumption & HDS

Hydrogen consumption (NL/L) vs. TOS



- Hydrogen consumption for SHU catalysts close to experimental error.
- HDS catalysts could be clearly differentiated in terms of H<sub>2</sub> consumption which correlated with catalyst activity


c(S) vs TOS



- Top performing catalysts reduced total sulfur below 20 ppm at the highest temperature of the screening program
- For the final test all catalysts were brought close to 20 ppm total sulfur concentration via individual reactor temperature control
- SHU reactors with slightly higher S in product concentration due to evaporation of light ends (quantified via GC)

# Results


## Sulfur speciation



- Sulfur speciation achieved GC-SCD analysis and sample derivatization to remove Mercaptanes and Sulfides;
  - The derivatization process involves a two-step selective oxidation, followed by sample work-up for analysis.
- Bottom purple block for H<sub>2</sub>S still dissolved in sample



# Summary



- **HDS activity & Octane retention**
  - Benchmarking OctaMax for HDS activity and octane retention
  - Deep HDS with minimal olefin saturation
- **Operational efficiency**
  - Dedicated pilot plant testing under realistic conditions
- **Sulfur speciation performed in liquid products**
  - Product Distribution: PIONA analysis of both vapor and liquid products
- **Advanced analytics:**
  - Validation of S conversion, yield, mass balance and octane retention



hte

# Evonik in figures

## FINANCES

**€15.2**

**billion**

in sales generated by our company in the 2024 fiscal year.

## BUSINESS

**9,200**

**products**

in our portfolio – from ABIL® to NANOPOX® and ZETASPERSE®.

**104**

**production locations**

ensure close proximity to customers and markets, whether in North America, South America, Europe or Asia.

## INNOVATION

**21,400**

**patents**

stand for Evonik's innovative spirit. Our first patent dates back to 1882/83.

**€459**

**million**

invested in our company's research and development activities in 2024.


## PEOPLE

**32,000**

**employees**

represent 110 nations. Plenty of potential to develop tailor-made solutions for every market in the world.


# Hydroprocessing Catalyst: Products & Services






# OctaMax: Typical Gasoline Pool Composition

% Blend Stocks of Gasoline Pool




Distribution of Sulfur in Blend Stocks




# OctaMax: Sulfur and Olefin Distribution in FCC-Naphtha

Distribution of Sulfur




Distribution of Olefins



The goal is to minimize olefins saturation while deeply desulfurizing

# OctaMax: Typical Process Configuration of FCC Naphtha -HDS Units



## SHU

- Always a NiMo catalyst
- Also known as Di-Olefin Reactor

## HDS

- Always a CoMo catalyst for HDS (and Ni catalyst for mercaptan recombination)
- Hydrodesulfurization. Reduce sulfur levels to allow gasoline pool blend to be <10wppm

Uniquely selected catalyst  
regenerated at optimal  
conditions for FCC  
naphtha HDS units

.....


OctaMax™






# HDS vs. Temperature

- Competitor A, B, and C catalysts are the most commonly utilized FCC naphtha catalysts from three different manufacturers
- **OctaMax catalyst demonstrates significantly increased desulfurization activity**




# Selectivity

- OctaMax demonstrates equivalent octane selectivity to fresh competitive catalysts.



# OctaMax Candidates

Total Product Sulfur VS Temperature



Not all regenerated CoMo catalysts are equivalent for FCC naphtha desulfurization

All CoMo catalysts in the graph met Evonik's standard internal specifications for regeneration for reuse

CoMo A does not meet OctaMax specifications as it has significantly reduced desulfurization activity

# Testing confirmed performance projections in our Technical Model

## Able to Predict Typical Parameters

Able to predict typical parameters such as WABT requirements, hydrogen consumption, cycle length projections, and generate loading diagrams

Typical  
Parameters



## Ability to Predict RON

Uniquely for this application, we also have the ability to predict RON loss

- Can model FCC naphtha post-treaters with or without an interstage splitter


Predict  
RON



## Comparing Feed Slate Cases

With this information we can help our customers optimize run strategies by comparing different feed slate cases

Compare



Guarantee



## Guarantee

All modeling is backed by a performance guarantees consistent with industry peers



# Why this matters

## OctaMax

### Enhanced HDS activity

#### Advantages:

- Can treat more difficult feedstocks
- Increase feed rates
- Increase catalyst cycle lengths



### Equivalent selectivity toward octane retention

#### Advantages:

- Flexibility to overtreat when economics are favorable



### Lower fill cost

#### Advantages:

- Compared to fresh, OctaMax is a more cost effective catalyst



### Sustainability benefits

#### Advantages:

- Reduced CO<sub>2e</sub> compared to manufacturing fresh catalyst



# Why Testing with hte is valuable

1

Reliable,  
independent testing  
institution

2

Ability to test a large  
number of catalyst  
systems simultaneously

3

Experience and assets  
to accurately test many  
different applications

4

Willingness to customize test  
set-up and analytical  
methods to precisely fit  
customers needs

5

Flexibility to alter plans  
mid-test to best capture the  
intended data





**EVONIK**

Leading Beyond Chemistry

&

hte